

EXHIBIT F

(12) **United States Patent**
Nicolson et al.

(10) **Patent No.: US 6,951,894 B1**
(45) **Date of Patent: *Oct. 4, 2005**

(54) **EXTENDED WEAR OPHTHALMIC LENS**

(75) Inventors: **Paul Clement Nicolson**, Dunwoody, GA (US); **Richard Carlton Baron**, Alpharetta, GA (US); **Peter Chabreck**, Basel (CH); **John Court**, Ultimo (AU); **Angelika Domschke**, Lörrach (DE); **Hans Jörg Griesser**, The Patch (AU); **Arthur Ho**, Randwick (AU); **Jens Höpken**, Lörrach (CH); **Bronwyn Glenice Laycock**, Heidelberg Heights (AU); **Qin Liu**, Duluth, GA (US); **Dieter Lohmann**, Munchestein (CH); **Gordon Francis Meijs**, Murrumbena (AU); **Eric Papaspiliotopoulos**, Paddington (AU); **Judy Smith Riffle**, Blacksburg, VA (US); **Klaus Schindhelm**, Cherrybrook (AU); **Deborah Sweeney**, Roseville (AU); **Wilson Leonard Terry, Jr.**, Alpharetta, GA (US); **Jürgen Vogt**, Fribourg (CH); **Lynn Cook Winterton**, Alpharetta, GA (US)

(73) Assignee: **CIBA Vision Corporation**, Duluth, GA (US)

(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 0 days.

This patent is subject to a terminal disclaimer.

(21) Appl. No.: **09/640,526**

(22) Filed: **Aug. 17, 2000**

Related U.S. Application Data

(63) Continuation of application No. 09/262,542, filed on Mar. 4, 1999, now abandoned, which is a continuation of application No. 09/108,714, filed on Jul. 1, 1998, now Pat. No. 5,965,631, which is a division of application No. 08/682,452, filed on Jul. 17, 1996, now Pat. No. 5,849,811, which is a division of application No. 08/569,816, filed on Dec. 8, 1995, now Pat. No. 5,760,100, which is a continuation-in-part of application No. 08/301,166, filed on Sep. 6, 1994, now abandoned.

(30) **Foreign Application Priority Data**

Apr. 4, 1995 (DE) 95810221
May 19, 1995 (CH) 1496/95

(51) **Int. Cl.** **C08K 3/00**

(52) **U.S. Cl.** **523/107; 523/106; 523/108**

(58) **Field of Search** **523/106, 107**

(56) **References Cited**

U.S. PATENT DOCUMENTS

2,718,516 A 9/1955 Bortnick
3,284,406 A 11/1966 Nelson
3,518,324 A 6/1970 Polmanteer
3,700,573 A 10/1972 Laizler
3,708,225 A 1/1973 Misch et al.

3,808,178 A 4/1974 Gaylord
3,916,033 A 10/1975 Merrill
3,935,342 A 1/1976 Lim
3,959,102 A 5/1976 Wajs et al.
3,959,105 A 5/1976 Feneberg et al.
3,996,187 A 12/1976 Travnick
4,008,198 A 2/1977 Krohberger et al.
4,062,627 A 12/1977 Wajs et al.
4,095,878 A 6/1978 Fanti
4,099,859 A 7/1978 Merrill
4,112,207 A 9/1978 Jones
4,114,993 A 9/1978 Travnick
4,120,570 A 10/1978 Gaylord
4,128,318 A 12/1978 Sieglaff et al.
4,136,250 A 1/1979 Mueller et al.
4,143,949 A 3/1979 Chen
4,152,508 A 5/1979 Ellis et al.
4,153,641 A 5/1979 Deichert et al.
4,156,066 A 5/1979 Gould
4,156,067 A 5/1979 Gould
4,166,255 A 8/1979 Graham
4,169,119 A 9/1979 Covington
4,182,822 A 1/1980 Chang
4,186,026 A 1/1980 Rotenberg et al.
4,189,546 A 2/1980 Deichert et al.
4,195,030 A 3/1980 Deichert et al.
4,197,266 A 4/1980 Clark et al.
4,198,131 A 4/1980 Birdsall et al.
4,208,362 A 6/1980 Deichert et al.
4,208,506 A 6/1980 Deichert et al.

(Continued)

FOREIGN PATENT DOCUMENTS

EP 0114894 A1 7/1983
EP 0108886 A2 9/1983

(Continued)

OTHER PUBLICATIONS

Plaintiff's Ciba Vision's Memorandum in Opposition to Defendant's Motion for Summary Judgment, May 31, 2002 with Exhibits 1-7.

Bausch & Lomb's Expert Reports, Jean M. Frechet, Jun. 24, 2002; Dr. Erwin A. Vogler, Jun. 24, 2002; Preliminary Report of William J. Benjamin, Jun. 24, 2002; Preliminary Report of Dr. William Daly, Jun. 24, 2002 and Nonpublicly Available Exhibits to Daly Report, Jun. 24, 2002.

(Continued)

Primary Examiner—Edward J. Cain

(74) *Attorney, Agent, or Firm*—Banner & Witcoff, Ltd.

(57) **ABSTRACT**

An ophthalmic lens suited for extended-wear periods of at least one day on the eye without a clinically significant amount of corneal swelling and without substantial wearer discomfort. The lens has a balance of oxygen permeability and ion or water permeability, with the ion or water permeability being sufficient to provide good on-eye movement, such that a good tear exchange occurs between the lens and the eye. A preferred lens is a copolymerization product of a oxypm macromer and an ionopm monomer. The invention encompasses extended wear contact lenses, which include a core having oxygen transmission and ion transmission pathways extending from the inner surface to the outer surface.

84 Claims, No Drawings

US 6,951,894 B1

Page 2

U.S. PATENT DOCUMENTS									
				4,687,816	A	*	8/1987	Lin et al.	525/279
				4,696,974	A		9/1987	Sulc et al.	
4,214,014	A	*	7/1980	Hofer					427/40
4,217,038	A		8/1980	Letter et al.			10/1987	Wingler et al.	
4,225,631	A		9/1980	Berger et al.			12/1987	Harvey, III	
4,228,269	A		10/1980	Loshaeck et al.			2/1988	Yamamoto et al.	
4,229,273	A		10/1980	Wajs			3/1988	Galin	
4,242,483	A		12/1980	Novicky			3/1988	Bawa et al.	
4,245,069	A		1/1981	Covington			4/1988	Bruns et al.	
4,254,248	A		3/1981	Friends et al.			4/1988	Falcetta et al.	
4,259,467	A		3/1981	Keogh et al.			4/1988	Su et al.	
4,260,725	A		4/1981	Keogh et al.			5/1988	Mizutani et al.	
4,261,875	A		4/1981	LeBoeuf			6/1988	Froix	
4,276,402	A		6/1981	Chromeczek et al.			9/1988	Ratkowski	
4,277,595	A		7/1981	Deichert et al.			10/1988	Su et al.	
4,291,953	A		9/1981	Covington			10/1988	Deichert	
4,294,974	A		10/1981	LeBoeuf			12/1988	Su et al.	
4,303,772	A		12/1981	Novicky			2/1989	Dunks et al.	
4,312,575	A		1/1982	Peyman et al.			2/1989	Goldgerg et al.	
4,322,517	A		3/1982	Deubzer et al.			3/1989	Friends et al.	
4,327,203	A		4/1982	Deichert et al.			4/1989	Rice et al.	
4,332,922	A		6/1982	Kossmehl et al.			5/1989	Eillis	
4,341,889	A		7/1982	Deichert et al.			5/1989	Stoyan	
4,355,135	A		10/1982	January			5/1989	Kunzler et al.	
4,355,147	A		10/1982	Deichert et al.			6/1989	Mueller et al.	
4,359,558	A		11/1982	Gould et al.			6/1989	Sweet et al.	
4,365,050	A		12/1982	Ivani			7/1989	Dillion	
4,395,496	A		7/1983	Wittmann et al.			8/1989	Schafer et al.	
RE31,406	E		10/1983	Gaylord			8/1989	Su et al.	
4,408,023	A		10/1983	Gould et al.			8/1989	Dillion	
4,410,674	A		10/1983	Ivani			9/1989	Geyer et al.	
4,413,104	A		11/1983	Deubzer et al.			10/1989	Froix	
4,423,195	A		12/1983	Covington			1/1990	Moreau et al.	
4,424,305	A		1/1984	Gould et al.			3/1990	Bambury et al.	
4,436,887	A		3/1984	Chromeczek et al.			4/1990	Schafer et al.	
4,439,583	A		3/1984	Gould et al.			5/1990	Mueller et al.	
4,439,584	A		3/1984	Gould et al.			5/1990	Khan et al.	
4,439,585	A		3/1984	Gould et al.			7/1990	Leach et al.	
4,440,918	A		4/1984	Rice et al.			7/1990	Frances et al.	
4,447,562	A		5/1984	Ivani			7/1990	Deichert et al.	
4,454,295	A		6/1984	Wittmann et al.			7/1990	Markel et al.	
4,454,309	A		6/1984	Gould et al.			8/1990	Wallsten et al.	
4,463,149	A		7/1984	Ellis			8/1990	Novicky	
4,478,981	A		10/1984	Arkles			8/1990	Fleischmann et al.	
4,486,577	A		12/1984	Mueller et al.			9/1990	Toyoshima et al.	
4,487,905	A		12/1984	Mitchell			9/1990	Mueller	
4,495,361	A		1/1985	Friends et al.			10/1990	Goldbert et al.	
4,496,535	A		1/1985	Gould et al.			10/1990	Harisiades	
4,500,676	A		2/1985	Balazs et al.			11/1990	Janssen et al.	
4,527,293	A		7/1985	Eckstein et al.			12/1990	Culberson et al.	
4,528,301	A		7/1985	Upchurch			12/1990	Sugiyama	427/444
4,543,398	A		9/1985	Bany et al.			1/1991	Hahn et al.	
4,546,123	A		10/1985	Schäfer et al.			1/1991	Mueller et al.	
4,550,139	A		10/1985	Arkles			3/1991	Stoyan	
4,555,372	A		11/1985	Kunzler et al.			4/1991	Kunzler et al.	
4,563,565	A		1/1986	Kampfer et al.			4/1991	Lee et al.	
4,576,850	A		3/1986	Martens			4/1991	Mueller	
4,582,884	A		4/1986	Ratkowski			4/1991	Mueller	
4,582,885	A		4/1986	Barber			5/1991	Piskoti	
4,602,074	A		7/1986	Mizutaki et al.			5/1991	Spinelli	
4,605,712	A		8/1986	Mueller et al.			6/1991	Onozuka et al.	
4,616,045	A		10/1986	Upchurch			7/1991	Baron et al.	
4,625,007	A		11/1986	Ellis et al.			7/1991	Lai	525/100
4,626,292	A		12/1986	Sherman			8/1991	Molock et al.	
4,632,844	A		12/1986	Yanagihara et al.			10/1991	Pinchik	
4,632,968	A		12/1986	Yokota et al.			10/1991	Spinelli	
4,649,184	A		3/1987	Yoshikawa et al.			11/1991	Wu et al.	
4,649,185	A		3/1987	Takamizawa et al.			12/1991	Robertson et al.	
4,652,622	A		3/1987	Friends et al.			12/1991	Robertson et al.	
4,659,777	A		4/1987	Riffle et al.			12/1991	Bambury et al.	
4,661,573	A		4/1987	Ratkowski et al.			12/1991	Broderick	
4,663,409	A		5/1987	Friends et al.			12/1991	Nordan	
4,664,657	A		5/1987	Williamitis et al.			12/1991	Schwabe et al.	
4,665,145	A		5/1987	Yokota et al.			1/1992	Druskis et al.	
4,666,249	A		5/1987	Bauman et al.			1/1992	Kindt-Larsen	
4,668,558	A		5/1987	Barber			1/1992	Kamel et al.	
							1/1992	Stoyan	

US 6,951,894 B1

Page 3

5,091,204	A	2/1992	Ratner et al.	
5,094,876	A	3/1992	Goldberg et al.	
5,098,618	A	3/1992	Zelez	
5,100,689	A	3/1992	Goldberg et al.	
5,104,213	A	4/1992	Wolfson	
5,106,930	A	4/1992	Gupta	
5,115,056	A	5/1992	Mueller et al.	
5,116,369	A	5/1992	Kushibiki et al.	
5,128,408	A	7/1992	Tanaka et al.	
5,128,434	A	7/1992	Lai	
5,135,297	A	8/1992	Valint, Jr.	
5,141,748	A	8/1992	Rizzo	
5,147,396	A	9/1992	Kageyama et al.	
5,154,861	A	10/1992	McBrierty et al.	
5,157,093	A	10/1992	Harisiades et al.	
5,158,573	A	10/1992	Berg	
5,158,717	A	* 10/1992	Lai 264/1.1	
5,160,597	A	11/1992	Colapicchioni et al.	
5,162,396	A	11/1992	Hilty et al.	
5,162,469	A	11/1992	Chen	
5,166,298	A	11/1992	Friedmann et al.	
5,171,607	A	12/1992	Cumbo	
5,171,809	A	12/1992	Hilty et al.	
5,177,165	A	1/1993	Valint, Jr. et al.	
5,177,167	A	1/1993	Tone et al.	
5,196,493	A	3/1993	Gruber et al.	
5,219,965	A	* 6/1993	Valint 526/245	
5,236,969	A	8/1993	Kunzler et al.	
5,238,613	A	8/1993	Anderson	
5,260,000	A	11/1993	Nandu et al.	
5,260,001	A	11/1993	Nandu et al.	
5,264,161	A	11/1993	Druskis et al.	
5,270,418	A	12/1993	Kunzler et al.	
5,271,875	A	12/1993	Appelton et al.	
5,274,008	A	12/1993	Lai	
5,298,533	A	3/1994	Nandu et al.	
5,310,779	A	5/1994	Lai	
5,314,960	A	5/1994	Spinelli et al.	
5,321,108	A	6/1994	Kunzler et al.	
5,334,681	A	* 8/1994	Mueller 526/243	
5,336,797	A	8/1994	McGee et al.	
5,346,946	A	9/1994	Yokoyama et al.	
5,352,714	A	10/1994	Lai et al.	
5,357,013	A	10/1994	Bambury et al.	
5,358,995	A	10/1994	Lai et al.	
5,364,918	A	11/1994	Valint, Jr. et al.	
5,371,147	A	12/1994	Spinelli et al.	
5,374,662	A	12/1994	Lai et al.	
5,378,412	A	1/1995	Smith et al.	
5,387,632	A	2/1995	Lai et al.	
5,387,662	A	2/1995	Kunzler et al.	
5,387,663	A	2/1995	McGee et al.	
5,391,589	A	* 2/1995	Kiguchi 523/106	
5,420,324	A	5/1995	Lai et al.	
5,435,943	A	7/1995	Adams et al.	
5,449,729	A	9/1995	Lai et al.	
5,451,617	A	9/1995	Lai et al.	
5,451,651	A	9/1995	Lai	
5,453,467	A	9/1995	Bamford et al.	
5,456,864	A	10/1995	Wickes et al.	
5,466,147	A	11/1995	Appelton et al.	
5,486,579	A	1/1996	Lai et al.	
5,496,871	A	3/1996	Lai et al.	
5,512,205	A	4/1996	Lai	
5,525,691	A	6/1996	Valint, Jr. et al.	
5,534,605	A	7/1996	Bambury et al.	
5,539,016	A	7/1996	Kunzler et al.	
5,563,184	A	10/1996	McGee et al.	
5,610,252	A	3/1997	Bambury et al.	
5,629,360	A	5/1997	Askari et al.	
5,726,733	A	3/1998	Lai et al.	
5,849,811	A	* 12/1998	Nicolson et al. 523/106	
5,760,100	A	11/2000	Nicolson et al.	
5,776,999	A	11/2000	Nicolson et al.	
5,789,461	A	11/2000	Nicolson et al.	
5,849,811	A	11/2000	Nicolson et al.	

FOREIGN PATENT DOCUMENTS

EP	0184729	B1	11/1985
EP	0277771	A3	1/1988
EP	0 295 947	B1	6/1988
EP	0295947	A3	6/1988
EP	0170141	B1	7/1988
EP	0306756	A2	8/1988
EP	0306756	B1	8/1988
EP	0330616	A1	2/1989
EP	0395583	A2	4/1990
EP	90810308		4/1990
EP	0425436	A3	10/1990
EP	0461270	A1	12/1990
EP	0489185	A1	12/1990
EP	0439429	A2	1/1991
EP	0 443 005		4/1991
EP	0480238	A1	9/1991
EP	0497204	A2	1/1992
EP	0528664	A1	2/1993
EP	0584764	A1	8/1993
EP	0614921	A2	3/1994
EP	0620455	A2	4/1994
EP	0687550	A2	6/1994
EP	0643083	A1	9/1994
EP	0686491	A2	6/1995
EP	0735097	A1	3/1996
EP	0781777	A1	12/1996
EP	0782016	A2	12/1996
JP	6[1994]-122779		5/1994
WO	WO90/09013		8/1990
WO	WO92/07013		4/1992
WO	WO 94/15980		7/1994
WO	WO94/15980		7/1994
WO	WO95/24187		9/1994
WO	WO95/17689		6/1995
WO	WO95/20476		8/1995

OTHER PUBLICATIONS

Bausch & Lomb's Incorporated's Reply Brief In Support of Its Motion For Summary Judgment, Jun. 20, 2002 with Exhibits A-I.

Defendant Bausch & Lomb Incorporated's Motion for Summary Judgment of Invalidity Pursuant to 35 U.S.C. § 112, ¶ 1., May 7, 2002 with Exhibits A-M.

Defendant Bausch & Lomb Incorporated's Memorandum in Support of its Motion for Summary Judgment of Invalidity Pursuant to 35 U.S.C. § 112, ¶ 1., May 7, 2002.

Trial vol. 1, *Wesley Jessen Corporation v. Bausch & Lomb, Inc.*, May 29, 2002; Trial vol. 2, May 30, 2002; Trial vol. 3, May 31, 2002; Trial vol. 4, Jun. 3, 2002; and Trial vol. 5, Jun. 4, 2002.

Expert Report of Robert L. Sheridan with Exhibits A, B and C.

Carmaels & Ransford, documents relating to opposition—Jun. 10, 2002.

Ciba's Surreponse in Opposition to B&L's Motion for Summary Judgment on Best Mode Grounds, Jul. 3, 2002 with Exhibit A, Deposition of Lynn Cook Winterton and Exhibit B, Deposition of Richard Baron.

United States Court of Appeals, Reply Brief of Appellant Bausch & Lomb Incorporated, dated Oct. 16, 2002.

United States District Court for the District of Delaware, Decision dated Jun. 26, 2002 on *Wesley Jessen Corporation v. Bausch & Lomb, Inc.*

United States Court of Appeals for the Federal Circuit, Brief of Appellant Bausch & Lomb Incorporated, dated Aug. 22, 2002.

United States Court of Appeals for the Federal Circuit, Brief of Appellee Wesley Jessen Corporation, dated Oct. 4, 2002.

US 6,951,894 B1

Page 4

Plaintiff's Ciba Vision's Memorandum in Opposition to Defendant's Motion for Summary Judgment, May 31, 2002 with Exhibits 1-7. (Note: This Information may be Subject to Protective Order—not Open to Public, Only by Examiner or Other Authorized Patent and Trademark Office Employee).

Bausch & Lomb's Expert Reports, Jean M. Frechet, Jun. 24, 2002; Dr. Erwin A. Vogler, Jun. 24, 2002; Preliminary Report of William J. Benjamin, Jun. 24, 2002. (Note: This Information may be Subject to Protective Order—not Open to Public, only by Examiner or other Authorized Patent and Trademark Office Employee).

Bausch & Lomb's Incorporated's Reply Brief In Support of Its Motion For Summary Judgment, Jun. 20, 2002 with Exhibits A-I. (Note: This Information may be Subject to Protective Order—Not Open to Public, Only by Examiner or Other Authorized Patent and Trademark Office Employee). Defendant Bausch & Lomb Incorporated's Motion for Summary Judgment of Invalidity Pursuant to 35 U.S.C. § 112, ¶ 1., May 7, 2002 with Exhibits A-M. (Note: This Information may be Subject to Protective Order—Not Open to Public, Only by Examiner or Other Authorized Patent and Trademark Office Employee).

Defendant Bausch & Lomb Incorporated's Memorandum in Support of its Motion for Summary Judgment of Invalidity Pursuant to 35 U.S.C. § 112, ¶ 1., May 7, 2002. Note: This Information may be Subject to Protective Order—Not Open to Public, Only by Examiner or Other Authorized Patent and Trademark Office Employee).

Trial Volume 1, Wesley Jessen Corporation v. Bausch & Lomb, Inc., May 29, 2002; Trial Volume 2, May 30, 2002; Trial Volume 3, May 31, 2002; Trial Volume 4, Jun. 3, 2002; and Trial Volume 5, Jun. 4, 2002.

Expert Report of Robert L. Sheridan with Exhibits A, B and C. (Note: This Information may be Subject to Protective Order—not Open to Public, Only by Examiner or Other Authorized Patent and Trademark Office Employee).

Carmaels & Ransford, documents relating to opposition—Jun. 10, 2002.

Ciba's Surreponse in Opposition to B&L's Motion for Summary Judgment on Best Mode Grounds, Jul. 3, 2002 with Exhibit A, Deposition of Lynn Cook Winterton and Exhibit B, Deposition of Richard Baron. (Note: This Information may be Subject to Protective Order—not Open to Public, Only by Examiner or Other Authorized Patent and Trademark Office Employee).

Declaration of Yu-Chin Lai—Dec. 18, 2001.

Bausch & Lomb's Response to Ciba's Surreply Brief—Jan. 15, 2002.

Designated "Highlighted" Pages of Deposition of Yu-Chin-Lai, Nov. 17, 1999 (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

R & D Clinical Report, with "Highlighted" Portions, Feb. 7, 1994 (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

Designated "Highlighted" Pages of Deposition of Paul Nicolson, Nov. 10, (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

Distribution with "Highlighted" Pages (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

Bausch & Lomb Incorporated's Opposition to CIBA Corporation's Motion to Permit Certain Attorney's to have access to Information Under the Stipulated Protective Order (with Exhibits A-J) (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

Defendant Bausch & Lomb Incorporated's Opposition to CIBA Corporation's Motion to Use Evidence from this case at the Reexamination Proceedings (with Exhibits A-J) (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

An article in *Macromolecular Engineering* (1995) entitled "*Transparent Multiphasic Oxygen Permeable Hydrogels based on Silaxanic Statistical Copolymers*" by Robert et al, which was publicly available in Aug. 1995 in the United States (Exhibit G).

An article in *Trends in Polymer Science* entitled "*Silicone Hydrogels for Contact Lens Application*" by Jay F. Kunzler which was publicly available in Feb. 1996 in the United Kingdom (Exhibit H).

An article in the *Contact Lens Association of Ophthalmologists Journal* ("CLAO"), vol. 13, No. 1 (Jan. 1987) entitled "*Gas Permeable Hard Contact Lens Extended Wear: Ocular and Visual Responses to a 6-Month Period of Wear*" by Kenneth A. Polse et al., which was publicly available in Jan. 1987 in the United States of America (Exhibit I).

An article in the *Contact Lens Association of Ophthalmologists Journal* ("CLAO"), vol. 15, No. 2 (Apr. 1989) entitled "*Paraperm® EW Lens for Extended Wear*" by James E. Key et al., which was publicly available in Apr. 1989 in the United States of America (Exhibit J).

The publication of an article in *Biomaterials*, vol. 4 (Oct. 1983) entitled "*Gas-to-Liquid permeation in silicon-containing crosslinked, glassy copolymers of methyl methacrylate*" by Wu-Huang, Michael Yang, et al. which was publicly available in Oct. 1983 in the United States of America (Exhibit K).

The publication of an article in *Macromolecular Symposia*, vol. 98 entitled "*Recent Advances in the Design of Polymers for Contact Lenses*" by Gary D. Friends, et al., which was publicly available in Jul. 1995 in the United States of America (Exhibit L).

An article in *The British Polymer Journal* (Dec. 1976) entitled "*Polymers in Contact Lens Applications VI, The 'Dissolved' Oxygen Permeability of Hydrogels and the Design of Materials for use in Continuous-Wear Lenses*" by Chiong O. Ng et al., which was publicly available in Dec. 1976 in the United Kingdom (Exhibit M).

An article in the *Journal of Japan Contact Lens Society*, vol. 21, No. 5 (1979) entitled "*Physical Properties of High Water Combined Lens Materials, Compared with Rabbit's Cornea*" by Yuji Kosaka et al., which was publicly available in 1979 in Japan (Exhibit N).

The Medical Device Approval Letter, vol. IV, No. 1, which was publicly available on line through a number of databases including DIALOG and LEXIS in Jan. 1995 (Exhibit O).

US 6,951,894 B1

Page 5

The Medical Devices, The Diagnostics & Instrumentation Report—“*The Grey Sheet*”, vol. 21, No. 4, dated Jan. 23, 1995 which was publicly available on line through a number of databases including DIALOG and LEXIS in about Jan. 1995 (Exhibit P).

The B&L US 510(k) application filed with the U.S. Food and Drug Administration (“USFDA”) for the Premier 90 (Balafilicon A) contact lens which application was publicly available from the USFDA from Jan. 1995 (Exhibit Q).

CIBA Vision Corporation, a Delaware corporation vs. Bausch & Lomb Incorporated, a New York Corporation—Counter Claims—United States District Court for the Northern District of Georgia (Gainesville)—pp. 22–28 with Exhibit F—Apr. 19, 1999.

CIBA Vision Corporation, a Delaware corporation vs. Bausch & Lomb Incorporated, a New York Corporation—Memorandum in Opposition to CIBA’s Motion for a Temporary Restraining Order and a Preliminary Injunction—United States District Court for the Northern District of Georgia (Gainesville) pp. 2–17 with Exhibits A–F and J. In the Federal Court of Australia Victoria District Registry—Defence and Cross-Claim—No. V81 of 2000—Jun. 2, 2000.

Releasable 510(k) No. K972454—Aug. 8, 1997—Lenses, Soft Contact, Daily Wear.

Releasable 510(k) No. K944895—Dec. 8, 1994—Contact Lenses.

Letter from Mr. Edward W. Remus to Michael O. Sutton, Esq.—Apr. 4, 2000—RE: *CIBA Vision Corporation v. Bausch & Lomb Incorporated*.

Letter from Mr. Edward W. Remus of McAndrews, Held & Malloy, Ltd. to Kenneth L. Cage, Esq. at McDermott, Will & Emery—RE: *CIBA Vision Corporation v. Bausch & Lomb Incorporated*—Jun. 5, 2000.

Letter from Richard T. McCaulley Jr. of McAndrews, Held & Malloy, Ltd., to Michael O. Sutton at Sidley & Austin—RE: *CIBA Vision Corp v. Bausch & Lomb Inc.* Civil Action No. 99–0034–WCO—Jun. 23, 2000.

Letter from Kenneth L. Cage of McDermott, Will & Emery to Edward W. Remus at AcAndrews, Held & Molloy, LTD.—RE: *CIBA Vision Corporation v. Bausch & Lomb Incorporated*—May 25, 2000.

Letter from Kenneth L. Cage of McDermott, Will & Emery to Colleague—RE: *CIBA Vision Corporation v. Bausch and Lomb, Incorporated* U.S. Patent Reexamination Proceedings U.S. Patent Nos. 5,760,100; 5,849,811; 5,789,461; and 5,776,999.

Letter from Carol Kirby of the University of New South Wales to Mr. R. Hearhard at CIBA Vision Corporation—Jul. 7, 2000—RE: Rule 56 Inquiry—UNSW CCLRU.

Letter from Michael O. Sutton of Sidley & Austin to Kenneth L. Cage of McDermott, Will & Emery and to Richard Gearhart at CIBA Vision Corp.—Jul. 12, 2000—RE: *CIBA Vision Corporation v. Bausch & Lomb Incorporated* (16008/0010).

Letter from Michael O. Sutton to Edward W. Remus—Dec. 10, 1999.

Letter from Michael O. Sutton to Edward W. Remus—Dec. 14, 1999.

Letter from Edward W. Remus to Michael O. Sutton—Dec. 17, 1999.

Letter from Michael O. Sutton to Edward W. Remus—Dec. 22, 1999.

Letter from Edward W. Remus to Michael O. Sutton—Jan. 3, 2000.

Letter from Michael O. Sutton to Edward W. Remus—Jan. 4, 2000.

Letter from Edward W. Remus to Michael O. Sutton—Jan. 6, 2000.

Letter from Michael O. Sutton to Edward W. Remus—Jan. 12, 2000.

Letter from Russell Wheatley to Edward W. Remus—Jan. 31, 2000.

Letter from Edward W. Remus to Michael O. Sutton—Feb. 2, 2000.

Letter from Edward W. Remus to Michael O. Sutton—Feb. 7, 2000.

Letter from Michael O. Sutton to Edward W. Remus—Feb. 10, 2000.

Letter from Edward W. Remus to Michael O. Sutton—Feb. 18, 2000.

Letter from Michael O. Sutton to Edward W. Remus—Mar. 23, 2000.

Letter from Edward W. Remus to Michael O. Sutton—Mar. 30, 2000.

Letter from Michael O. Sutton to Edward W. Remus—Apr. 3, 2000.

Letter from Michael O. Sutton to Edward W. Remus—Jun. 19, 2000.

Letter from Edward W. Remus to Michael O. Sutton—Jun. 30, 2000.

Weissman et al. “Cancellation of the Boundary and Edge Effects by Choice of Lens Thickness during Oxygen Permeability Measurement of Contact Lens”, *Optometry & Vision Science*, vol. 66, No. 5, pp. 264–268, (1989).

Hamano et al., “A study of the Complications Induced by Conventional and Disposable Contact Lenses”, *The CLAO Journal*, vol. 20, pp. 103–110, (1994).

Schein et al., “The Impact of Overnight Wear on the Risk of Contact Lens–Associated Ulcerative Keratitis”, *Arch. Ophthalmol*, vol. 112, pp. 186–190, (1994).

Nilsson et al., “The Hospitalized Cases of Contact Lens Induced Keratitis in Sweden and their Relation to Lens Type and Wear Schedule; Results of a Three–Year Retrospective Study”, *The CLAO Journal*, vol. 20, pp. 97–101, (1994).

Goyal et al., “Corneal Ulcers in Patients with Disposable Extended–Wear Contact Lenses”, *Ann. Ophthalmol—Glaucoma*, vol. 26, pp. 194–199, (1994).

Schein et al., “The Relative Risk of Ulcerative Keratitis Among Users of Daily–Wear and Extended–Wear Soft Contact Lenses”, *The New England Journal of Medicine*, vol. 321, pp. 773–776, (1989).

Solomon et al., “Testing Hypotheses for Risk Factors for Contact Lens–Associate Infectious Keratitis in an Animal Model”, *The CLAO Journal*, vol. 20, pp. 109–113, (1994).

Poggio et al., “Complications and Symptons in Disposable Extended Wear Lenses Compared with Conventional Sofy Daily Wear and Soft Extended Wear Lenses”, *The CLAO Journal*, vol. 19, pp. 31–39, (1993).

Poggio et al., “The Incidence of Ulcerative Keratitis Among Users of Daily–Wear and Extended–Wear Soft Contact Lenses”, *The New England Journal of Medicine*, vol. 321, pp. 779–783, (1989).

McRae et al., “Corneal Ulcer and Adverse Reaction Rates in Premarket Contact Lens Studies”, *American Journal of Ophthalmology*, vol. 111, pp. 457–465, (1991).

US 6,951,894 B1

Page 6

Grant et al., "Extended Wear of Hydrogel Lenses", Problems in Optometry, vol. 2, pp. 599-622, (1990).

Brennan et al., "Extended Wear in Perspective", Optometry and Vision Science, vol. 74, 609-623, (1997).

Winterton et al., "Coulometric Method for Measuring Oxygen Flux and Dk of Contact Lenses and Lens Materials", vol. 14, No. 11, 1987, pp. 441-449, (1987).

Winterton et al., "Coulometrically Determined Oxygen Flux and Resultant Dk of Commercially Available Contact Lenses", vol. 15, pp. 117-123, (1988).

Holden et al., "The Dk Project: An Interlaboratory Comparison of Dk/L Measurements", Optometry and Vision Science, vol. 67, pp. 476-481, (1989).

Winterton et al., "Coulometric Method for Measuring Oxygen Flux and Dk of Contact Lenses and Lens Materials", The Cornea: Transactions of The World Congress on The Cornea III, Raven Press, Chapter 47, pp. 273-280, (1988).

Alvord et al., "Oxygen Permeability of a New Type of High Dk Soft Contact Lens Material", Optometry and Vision Science, vol. 75, pp. 30-36, (1998).

Holden et al., "Critical Oxygen Levels to Avoid Corneal Edema for Daily and Extended Wear Contact Lenses", Investigative Ophthalmology & Visual Science, vol. 25, pp. 2-8, (1983).

Domschke et al., "Morphology Requirements For On-Eye Mobility of Soft Oxygen Permeable Contact Lenses", Business and Technology Innovation, Ciba Vision.

Kirk Othmer Encyclopedia of Chemical Technology (chapter on contact lenses), pp. 720-741.

Edelhauser et al. "Oxygen, calcium, lactate, chloride, and sodium permeability of soft contact lenses" in Gasset et al. "Soft Contact Lenses", The C.V. Mosby Co., (1972) Saint Louis, pp. 37-43.

U.S. Prosecution History of U.S. Pat. No. 5,034,461.

U.S. Prosecution History of U.S. Pat. No. 5,451,617.

U.S. Prosecution History of U.S. Pat. No. 5,486,579.

Mandell, R., *Contact Lens Practice, Fourth Edition*, Springfield, IL, pp. 520-521 (1988).

Lai, Yu-Chin, "The Role of Bulky Polysiloxanylalkyl Methacrylates in Oxygen Permeable Hydrogel Materials," ACS Workshop-Hyphenated Methods In Polymer Spectroscopy, *Polymeric Materials Science and Engineering*, vol. 69, Fall Meeting, Chicago, IL, pp. 228-231 (1993).

Fatt, I., et al., "Measuring Oxygen Permeability of Gas Permeable Hard and Hydrogel Lenses and Flat Samples in Air," vol. 14, No. 10 (Oct. 1987) pp. 389-401 (Oct. 1987).

Lai, Yu, et al., "Surface Wettability Enhancement of Silicone Hydrogel Lenses by Processing with Polar Plastic Molds," *J. of Biomed. Mat. Res.*, vol. 35:3, pp. 349-356 (1997).

Lai, Yu, et al., "Synthesis and Characterization of α,ω -Bis(4-hydroxybutyl) Polydimethylsiloxane," *Polymer Preprints*, 35:2, American Chemical Society (Aug. 1994).

Memorandum in Support of CIBA's Motion to Use Evidence from this Case in the Reexamination Proceedings ("CIBA Motion") (with Exhibits 1-8) (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

Reply Brief in Support of CIBA VISION Corporation's Motion to use Evidence from this Case in the Reexamination Proceedings ("CIBA Reply") (with Exhibits 1-7) (Note: This Information is Subject to a Protective Order—It is not Open to the Public, but Only to the Examiner or Other Authorized Patent and Trademark Office Employee).

The Journal of Clinical Pharmacology and Therapeutics, vol. 58, No. 5, p. 604. Nov. 1995 edition, which contained the structural formula for Balafilcon A and which was publicly available in Nov. 1995 in the United States (Exhibit A).

The publication in Chemistry and Industry of an article by Kunzler entitled "*Contact Lens Materials*", Aug. 21, 1995, pp. 651-655 and which was publicly available from Oct. 19, 1995 at the University of New South Wales Library (Exhibit B).

The publication of an article in the Journal of Applied Polymer Science, vol. 55 (1995), pp. 611-619, by Kunzler et al. entitled "*Hydrogels Based on Hydrophilic Side-Chain Siloxanes*", and which was publicly available from Feb. 21, 1995 at the University of New South Wales Library (Exhibit C).

An article in World Plastics and Rubber Technology (1992), pp. 23-29, entitled "*Clear Vision for the Future*" by Lai et al. and which was publicly available from Nov. 16, 1992 at the State Library of New South Wales (Exhibit D).

An article in the Journal of Applied Polymer Science, vol. 56 (1995) entitled "*Role of Bulky Polysiloxanylalkyl Methacrylates in Oxygen-Permeable Hydrogen Materials*" by Lai et al. and which was publicly available in Apr. 1995 in the United States (Exhibit E).

An article in the Journal of Applied Polymer Science, vol. 56 (1995) entitled "*Novel Polyurethane-Silicone Hydrogels*" by Lai which was publicly available in Apr. 1995 in the United States (Exhibit F).

Hung, William, Preparation and Testing Silicone Containing Hydrogel Lenses as Described in U.S. Pat. No. 52160000, submitted to EP patent office regarding EP Patent No. 0819258, Apr. 1, 2004.

Breuer & Muller, Grounds of Appeal, in EP Patent application 96 906 116.5-2217, dated Apr. 1, 2004.

Fisher, Adrian John, Submissions of Opponent II in response to the Proprietors' Appeal, dated Aug. 24, 2004.

Winterton, Lynn C., et al., Coulometric Method for Measuring Oxygen Flux and Dk of Contact Lenses and Lens Materials, ICLC, vol. 14, No. 11, Nov. 1987, pp. 441-450.

Declaration of William J. Benjamin, Professor of Optometry and Physiological Optic University of Alabama at Birmingham, Aug. 24, 2004.

Package insert, Focus Night & Day Soft Contact Lenses for Extended Wear, print date: Sep. 2001.

Order of J. Merkel in Federal Court of Australia, Victoria District Registry, Novartis AG v. Bausch & Lomb (Australia) Pty. Ltd [2004] FCA 835, dated Aug. 24, 2004.

Criticism of the Experimental Report provided by Dr. William M. Hung, submitted to EP patent office regarding EP Patent No. 0819258, Aug. 24, 2004.

Heidelberg, Alfred Huthig Verlag; Praktikum Der Makromolekularen Organischen Chemie, 1996. (Translation of relevant section included).

US 6,951,894 B1

Page 7

EP Patent application No. 90810308.8; Answer to EPO concerning Office Action, dated Mar. 30, 1994. (Translation of relevant section included).

Ciba Vision Corporation et al. v. Bausch & Lomb Incorporated, Plaintiff CIBA Vision Corporation's Post-Trial Proposed Findings of Fact and Conclusions of Law, United States District Court for the Northern District of Georgia (Gainesville), Civil Action File No. 2:99-0034-RWS, May 7, 2004, pp. 1-197.

Ciba Vision Corporation et al. v. Bausch & Lomb Incorporated, Bausch & Lomb's Proposed Findings of Fact and Conclusions of Law, United States District Court for the Northern District of Georgia (Gainesville), Civil Action File No. 2:99-0034-RWS, May 5, 2004, pp. 1-242.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Consent Decree and Order of Dismissal with Prejudice, United States District Court for the Northern District of Georgia, Gainesville Division, Civil Action No. 2:99-0034-RWS, Jul. 23, 2004.

Novartis AG v. Bausch & Lomb (Australia) Pty Ltd [2004] FCA 835, Consent Order, Federal Court of Australia, Victoria District Registry, Aug. 24, 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 1-33, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 34-67, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 68-138, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 141-274, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 276-422, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 423-667, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 675-782, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 787 B, C & D, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States

District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 789-844, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 846-941, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 944-1015, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 1016-1081, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 1082-1132, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Trial Exhibits 1134-1188, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 1-35, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 36-93, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibit 94, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 95-183, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 184-332, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 353-567, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 601-619, Apr. 2004.

US 6,951,894 B1

Page 8

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibit 620, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Trial Exhibits 621–640, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, CIBA's Admitted Demonstrative Trial Exhibits 100–114, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Demonstrative Trial Exhibits 27–82, Apr. 2004.

CIBA Vision Corporation et al. v. Bausch & Lomb, Inc., Trial Before the Honorable Richard W. Story, United States District Court Judge for the Northern District of Georgia, Atlanta Division, Bausch & Lomb's Admitted Demonstrative Trial Exhibit 83, Apr. 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #01 & #02 **BITC** CibaVision Jun. 30, 1993, CV51 002700.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #03 & #04 **BITC** CibaVision Jun. 30, 1993, CV51 002701.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #05 & #06 **BITC** CibaVision Jun. 30, 1993, CV51 002702.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #07 & #08 **BITC** CibaVision Jun. 30, 1993, CV51 002703.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #09 & #10 **BITC** CibaVision Jun. 30, 1993, CV51 002704.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #11 & #12 **BITC** CibaVision Jun. 30, 1993, CV51 002705.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #13 & #14 **BITC** CibaVision Jun. 30, 1993, CV51 002707.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #15 & #16 **BITC** CibaVision Jun. 30, 1993, CV51 002709.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #17 & #18 **BITC** CibaVision Jun. 30, 1993, CV51 002710.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #19 & #20 **BITC** CibaVision Jun. 30, 1993, CV51 002711.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #21 & #22 **BITC** CibaVision Jun. 30, 1993, CV51 002712.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #23 & #24 **BITC** CibaVision Jun. 30, 1993, CV51 002713.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #25 & #26 **BITC** CibaVision Jun. 30, 1993, CV51 002708.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #27 & #28 **BITC** CibaVision Jun. 30, 1993, CV51 002706.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision "See 3" Tape #29 RP Golf "Awards–Tape" CibaVision Jun. 30, 1993, CV51 002714.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Tape, 1. See3 Z–Por/Dextran Jun. 10, 1995 Time—1:50 2. See3 Dispensing Study Jun. 10, 1995 Infiltrates—Mcov/Raus Time 2:00, CV51 002531.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Tape, Extended Wear Brian Holden Ph.D, Sep. 16, 1992, Erba, Italy, CV51 002536.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision/CRU/CECER Research Collaboration Meeting, Sep. 16–20, 1992, Erba, Italy Tape 1 of 4 Time: 6 Hours, CV51 002532.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision/CRU/CECER Research Collaboration Meeting, Sep. 16–20, 1992, Erba, Italy Tape 2 of 4 Time: 6 Hours, CV51 002533.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision/CRU/CECER Research Collaboration Meeting, Sep. 16–20, 1992, Erba, Italy Tape 3 of 4 Time: 6 Hours, CV51 002534.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Ciba Vision/CRU/CECER Research Collaboration Meeting, Sep. 16–20, 1992, Erba, Italy Tape 4 of 4 Time: 6 Hours, CV51 002535.

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 7, 2003, (pp. P–1–P–76).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 8, 2003, (pp. P–77–P–153).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd. Oct. 9, 2003, (pp. P–154–P–224).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 10, 2003, (pp. P–225–P–294).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 13, 2003, (pp. P–295–P–359).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 14, 2003, (pp. P–360–P–422).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 15, 2003, (pp. P–422–P–484).

US 6,951,894 B1

Page 9

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 16, 2003, (pp. P-485-P-563).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 17, 2003, (pp. P-564-P-644).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 20, 2003, (pp. P-645-P-679).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 21, 2003, (pp. P-679-P-750), (pp. P1-P5).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 22, 2003, (pp. P-751-P-822).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 23, 2003, (pp. P-823-895).

Spark and Cannon, Transcript of Proceedings, Merkel J., Federal Court of Australia, Hearing & Motion, Novartis AG and Another and Bausch and Lomb (Australia) Pty Ltd., Oct. 24, 2003, (pp. P-896-968).

Order ("construction of the patent claims")—U.S. Court for the Northern District of Georgia Gainesville Division, Civil Action 2:99-CV-0034-RWS, dated Mar. 14, 2003.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. I, pp. 1-195, Mar. 15, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. II, pp. 196-406, Mar. 16, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. III, pp. 407-596, Mar. 17, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. IV, pp. 599-809, Mar. 19, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. V, pp. 810-1039, Mar. 22, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. VI, pp. 1040-1251, Mar. 23, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. VII, pp. 1252-1423, Mar. 24, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. VIII, pp. 1424-1577, Mar. 26, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. IX, pp. 1578-1781, Mar. 29, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. X, pp. 1782-1999, Mar. 30, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XI, pp. 2000-2200, Mar. 31, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XII, pp. 2201-2404, Apr. 1, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XIII, pp. 2405-2548, Apr. 2, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XIV, pp. 2549-2769, Apr. 13, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XV, pp. 2770-2989, Apr. 14, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XVI, pp. 2990-3204, Apr. 19, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XVII, pp. 3205-3432, Apr. 20, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XVIII, pp. 3433-3626, Apr. 21, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XIX, pp. 3627-3865, Apr. 22, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XX, pp. 3866-4072, Apr. 23, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XXI, pp. 4073-4303, Apr. 26, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XXII, pp. 4304-4319, Apr. 27, 2004.

CIBA Vision Corporation v. Bausch & Lomb, Inc., Transcript of Trail Before the Honorable Robert W. Story, United State District Court Judge for the Northern District of Georgia, Atlanta Division, Vol. XXIII, pp. 4320-4456, Apr. 28, 2004.

US 6,951,894 B1

Page 10

Opposition of European Patent No. 0819258, Novartis AG, et al. by Johnson & Johnson Vision Care, Inc., Jun. 10, 2002.

Novartis AG., Bausch & Lomb Inc., Johnson & Johnson Vision Care, Inc., European Patent Office, Minutes of the Oral Proceedings before the Opposition Division, Sep. 17, 2003.

Johnson & Johnson Pacific Pty Limited v Novartis AG et al., In the Federal Court of Australia New South Wales District Registry, APPLICATION to revoke Australian Patent No. 704749, pp. 1–3, Sep. 11, 2003.

Johnson & Johnson Pacific Pty Limited v Novartis AG et al., In the Federal Court of Australia New South Wales District Registry, Particulars of Invalidity, pp. 1–10, Sep. 11, 2003.

Johnson & Johnson Pacific Pty Limited v Novartis AG et al., In the Federal Court of Australia New South Wales District Registry, Statement of Claim, pp. 1–2, Sep. 11, 2003.

Johnson & Johnson (New Zealand) Limited v. Novartis AG et al. In the High Court of New Zealand Wellington Registry, Respondents' Reply to the Applicant's Particulars of Objection as to Validity of New Zealand Letters Patent No. 304321, pp. 1–10, Jan. 23, 2004.

Johnson & Johnson (New Zealand) Limited v. Novartis AG et al. In the High Court of New Zealand Wellington Registry, Respondents' Reply to the Applicant's Particulars of Objection as to Validity of New Zealand Letters Patent No. 334977, pp. 1–9, Jan. 23, 2004.

Johnson & Johnson (New Zealand) Limited v. Novartis AG et al. In the High Court of New Zealand Wellington Registry, Notice of Application by Applicant for Directions Regarding the Proceeding, pp. 1–2, Oct. 28, 2003.

Johnson & Johnson (New Zealand) Limited v. Novartis AG et al. In the High Court of New Zealand Wellington Registry, Particulars of Objection as to Validity of New Zealand Letters Patent No. 334977, pp. 1–10, Oct. 28, 2003.

Johnson & Johnson (New Zealand) Limited v. Novartis AG et al. In the High Court of New Zealand Wellington Registry, Particulars of Objection as to Validity of New Zealand Letters Patent No. 304321, pp. 1–10, Oct. 28, 2003.

Johnson & Johnson (New Zealand) Limited v. Novartis AG et al. In the High Court of New Zealand Wellington Registry, Originating Application for Revocation of New Zealand Letters Patent Nos. 304321 and 334977, pp. 1–3, Oct. 28, 2003.

Bausch & Lomb's 1994 Annual Report p. 14.

George L. Grobe III, et al., Surface Characterization of an Oxidized Silicone Hydrogel, 1994, pp. 43–47.

Surfaces in BioMaterials '94, Symposium Notebook, Sep. 7–10, 1994, Scottsdale, Arizona.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Bausch & Lomb's Supplementary Submission Concerning the Morphology of the Purevision Product, No. V81 of 2000, pp. 1–7, Dec. 17, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Respondent/Cross-Claimants' Submissions in Opposition to Amendment, No. V81 of 2000, pp. 1–19, Aug. 22, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia New South Wales District Registry, Bausch & Lomb's Submissions in Answer, No. V81 of 2000, pp. 1–29, Dec. 5, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia New South Wales District Registry, Bausch & Lomb's Response to the Applicants' Supplementary Closing Submissions, No. V81 of 2000, pp. 1–22, Dec. 5, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Statement of Grounds in Support of Amendment, No. V81 of 2000, pp. 1–4, Jun. 6, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Respondent/Cross-Claimants' Grounds of Opposition to Amendment, No. V81 of 2000, pp. 1–8, Jul. 4, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Bausch & Lomb's Contention of Fact and Law, No. V81 of 2000, pp. 1–29, Sep. 29, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Applicants' Opening Submissions, No. V81 of 2000, pp. 1–16, Sep. 29, 2003 (Glossary of Ophthalmic Terms pp. 1–8).

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Applicants' Supplementary Closing Submissions, No. V81 of 2000, pp. 1–3, Nov. 28, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Applicants' Reply Submissions, No. V81 of 2000, pp. 1–12, Dec. 5, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Applicants' Closing Submissions, No. V81 of 2000, pp. 1–27, Nov. 19, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Applicants' Outline of Submissions for Amendment, No. V81 of 2000, pp. 1–5, Aug. 25, 2003.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Applicants' Further Supplementary Submissions, No. V81 of 2000, pp. 1–17, Undated.

Novartis AG et al. v. Bausch & Lomb (Australia) Pty. Ltd. et al., In the Federal Court of Australia Victoria District Registry, Minutes of Consent Order, No. V81 of 2000, pp. 1–12, Dec. 2003.

Johnson & Johnson Vision Care, Inc. v. Ciba Vision Corporation, Complaint for Damages, Declaratory and Injunctive (sic) Relief and Demand for Jury Trial, United States District Court Middle District of Florida Jacksonville Division, Case No. 3:03–CV800–J–25TEM, Sep. 17, 2003, pp. 1–9.

* cited by examiner

US 6,951,894 B1

1

EXTENDED WEAR OPHTHALMIC LENS

This application is a continuation of Ser. No. 09/262,542, filed Mar. 4, 1999 now abandoned, which is a continuation of Ser. No. 09/108,714, filed Jul. 1, 1998, now U.S. Pat. No. 5,965,631 which is a divisional of application Ser. No. 08/682,452, filed Jul. 17, 1996, now U.S. Pat. No. 5,849,811 which is a divisional of application Ser. No. 08/569,816, filed Dec. 8, 1995 now U.S. Pat. No. 5,760,100 which is a continuation-in-part of U.S. application Ser. No. 08/301,166, filed on Sep. 6, 1994 now abandoned. Priority is also claimed 119 for German Application No. 95810221.2 filed on Apr. 4, 1995 and Swiss Application No. 1496/95 filed on May 19, 1995.

BACKGROUND OF THE INVENTION**1. Field of the Invention**

This invention relates broadly to lenses and polymeric materials useful in optic and ophthalmic arts. More specifically, this invention relates to polymeric materials and treatment processes useful in the manufacture of contact lenses. Still more specifically, this invention relates to contact lenses useful as extended-wear contact lenses.

2. Description of the Related Art

A wide variety of research has been conducted in the field of biocompatible polymers. The definition of the term "biocompatible" depends on the particular application for which the polymer is designed. In the field of ophthalmic lenses, and in particular in the field of contact lenses, a biocompatible lens may be generally defined as one which will not substantially damage the surrounding ocular tissue and ocular fluid during the time period of contact. The phrase "ophthalmically compatible" more appropriately describes the biocompatibility requirements of ophthalmic lenses.

One ophthalmic compatibility requirement for contact lenses is that the lens must allow oxygen to reach the cornea in an amount which is sufficient for long-term corneal health. The contact lens must allow oxygen from the surrounding air to reach the cornea because the cornea does not receive oxygen from the blood supply like other tissue. If sufficient oxygen does not reach the cornea, corneal swelling occurs. Extended periods of oxygen deprivation causes the undesirable growth of blood vessels in the cornea. "Soft" contact lenses conform closely to the shape of the eye, so oxygen cannot easily circumvent the lens. Thus, soft contact lenses must allow oxygen to diffuse through the lens to reach the cornea.

Another ophthalmic compatibility requirement for soft contact lenses is that the lens must not strongly adhere to the eye. Clearly, the consumer must be able to easily remove the lens from the eye for disinfecting, cleaning, or disposal. However, the lens must also be able to move on the eye in order to encourage tear flow between the lens and the eye. Tear flow between the lens and eye allows for debris, such as foreign particulates or dead epithelial cells, to be swept from beneath the lens and, ultimately, out of the tear fluid. Thus, a contact lens must not adhere to the eye so strongly that adequate movement of the lens on the eye is inhibited.

While there exist rigid gas permeable ("RGP") contact lenses which have high oxygen permeability and which move on the eye, RGP lenses are typically quite uncomfortable for the consumer. Thus, soft contact lenses are preferred by many consumers because of comfort. Moreover, a contact lens which may be continuously worn for a period of a day or more (including wear during periods of sleeping) requires comfort levels which exclude RGP lenses as popular extended-wear candidates.

2

In order to balance the ophthalmic compatibility and consumer comfort requirements in designing a daily wear soft contact lens, polymers and copolymers of 2-hydroxyethyl methacrylate (HEMA) were developed. These hydrophilic polymers move well on the eye and provide sufficient oxygen permeability for daily wear. Certain soft contact lenses have been approved by the FDA for extended wear periods of up to about 6 nights of overnight wear and seven days of daily wear. However, the consumer cannot safely and comfortably wear these poly(HEMA) lenses for extended periods of seven days or more, because the oxygen permeability is insufficient. True extended wear (i.e., seven days or more) of these lenses may result, at a minimum, in corneal swelling and development of surface blood vessels in the cornea.

In order to improve oxygen permeability, polymers containing silicone groups were developed. A variety of siloxane-containing polymers have been disclosed as having high oxygen permeability. For example, see U.S. Pat. Nos. 3,228,741; 3,341,490; 3,996,187; and 3,996,189. However, polysiloxanes are typically highly lipophilic. The properties (e.g., lipophilicity, glass transition temperature, mechanical properties) of known polysiloxanes has resulted in contact lenses which adhere to the eye, inhibiting the necessary lens movement. In addition, polysiloxane lipophilicity promotes adhesion to the lens of lipids and proteins in the tear fluid, causing a haze which interferes with vision through the lens.

There have been attempts to blend the desirable hydrophilic properties of hydrophilic polymers, formed from monomers such as HEMA, with the desirable oxygen permeability of polymers formed from siloxane-containing monomers. For example, see U.S. Pat. Nos. 3,808,178; 4,136,250; and 5,070,169. However, prior attempts at producing a true extended wear contact lens have been unsuccessful, either because of the effect of the extended-wear lens on corneal health or because the lens would not move on the eye. Thus, there remains a need for an ophthalmically compatible, transparent polymeric material which is suited to extended periods of continuous contact with ocular tissue and tear fluid.

OBJECTS AND SUMMARY OF THE INVENTION

An object of the invention is to provide a material having a balance of oxygen permeability, ion permeability, on-eye movement and tear exchange, all of which are sufficient for corneal health and wearer comfort during extended periods of continuous wear.

Another object of the invention is to provide an ophthalmic lens capable of extended continuous wear periods of at least 24 hours without substantial adverse impact on ocular health or consumer comfort, and more preferably, to provide a lens capable of continuous wear of 4 to 30 days or more without substantial adverse impact on ocular health or consumer comfort.

A further object of the invention is to provide an ophthalmic lens capable of extended continuous wear periods of at least 24 hours without substantial corneal swelling or consumer discomfort, and more preferably, to provide a lens capable of continuous wear of 4, 7, 14 or 30 days or more without substantial corneal swelling or consumer discomfort.

Yet another object of the invention is to provide methods of forming an extended-wear ophthalmic lens.

Still a further object of the invention is to provide methods of testing and classifying ophthalmic lenses as candidates for true extended-wear.

US 6,951,894 B1

3

These and other objects of the invention are met by the various embodiments described herein.

One embodiment of the invention is an ophthalmic lens, suited to extended periods of wear in continuous, intimate contact with ocular tissue and tear fluid. The lens displays a balance of oxygen permeability and ion permeability sufficient to maintain good corneal health, adequate movement of the lens on the eye and wearer comfort during extended wear periods. The lens is formed by polymerization, preferably copolymerization, of (a) at least one oxyperm polymerizable material which is capable of polymerizing to form a polymer having a high oxygen permeability; and (b) at least one ionoperm polymerizable material which is capable of polymerizing to form a polymer having a high ion permeability. Preferably, the lens includes a core polymeric material and ophthalmically compatible surfaces. In a preferred embodiment, the surface is more hydrophilic and lipophobic than the core polymeric material.

Another embodiment of the invention is a method of forming an ophthalmic lens having high oxygen permeability and high ion permeability. The method includes the step of forming a core material, having an inner surface and an outer surface, such that at least one pathway for ion transport and at least one pathway for oxygen transport are present from the inner to the outer surface. In a preferred embodiment, the method includes treating the surface of the lens to render the surface more hydrophilic than the core.

A further embodiment of the invention is an ophthalmic lens comprising a polymeric material which has a high oxygen permeability and a high ion or water permeability, the polymeric material being formed from at least one polymerizable material including (a) at least one oxyperm segment and (b) at least one ionoperm segment. The lens displays a balance of oxygen permeability and ion permeability sufficient to maintain good corneal health, adequate movement of the lens on the eye and wearer comfort during extended wear periods.

Yet another embodiment of the invention is a method of using a contact lens having both an oxygen transmission pathway and an ion transmission pathway from inner to outer surface as an extended wear lens. The method includes (a) applying the lens to the ocular environment and (b) allowing the lens to remain in intimate contact with the ocular environment for a period of at least 24 hours without substantial adverse impact on corneal health or wearer comfort. A preferred method includes additional steps of (c) removing the lens from the ocular environment; (d) disinfecting the lens; (e) applying the lens to the ocular environment; and (f) allowing the lens to remain in intimate contact with the ocular environment for a period of at least an additional 24 hours. In a preferred embodiment, the lens is worn for a continuous period of at least seven days without substantial adverse impact on corneal health or wearer comfort.

OUTLINE OF DESCRIPTION OF THE PREFERRED EMBODIMENTS

I. DEFINITION OF TERMS

II. CORE POLYMER AND LENS

- A. Oxyperm polymerizable materials
- B. Ionoperm polymerizable materials
- C. Weight ratio of oxyperm to ionoperm polymerizable materials
- D. Morphology
- E. Bulk Water Content

4

F. Ion and Water Permeability

1. Ionoflux Ion Permeability Measurements
2. Ionoton Ion Permeability Measurements
3. Hydrodell Water Permeability Measurements

G. Oxygen Permeability and Transmissibility

H. Mechanical On-eye Movement Parameters

1. Tensile Modulus and Short Relaxation Time
2. Tangent Delta
3. Parameter Combinations

I. Examples of suitable materials

1. Material "A"
2. Material "B"
3. Material "C"
4. Material "D"

III. OPHTHALMICALLY COMPATIBLE SURFACES

IV. UTILITY

A. Ophthalmic lenses

B. Contact lenses

V. METHODS OF USE AS EXTENDED-WEAR LENSES

VI. METHODS OF MANUFACTURE OF LENSES

DESCRIPTION OF THE PREFERRED EMBODIMENTS

One embodiment of the present invention is an ophthalmically compatible, transparent lens suited to extended periods of continuous contact with ocular tissue and tear fluids. A particularly preferred embodiment of the invention is an extended-wear vision correction lens suited for safe and comfortable long term wear without removal. In order to properly describe the invention and to delineate the metes and bounds of the claims, a set of basic terms will be defined at the outset.

I. DEFINITION OF TERMS

An "ophthalmic lens", as used herein, refers to lenses which are placed in intimate contact with the eye or tear fluid, such as contact lenses for vision correction (e.g., spherical, toric, bifocal), contact lenses for modification of eye color, ophthalmic drug delivery devices, ocular tissue protective devices (e.g., ophthalmic healing promoting lenses), and the like. A particularly preferred ophthalmic lens is an extended-wear contact lens, especially extended-wear contact lenses for vision correction.

A "polymerizable material which is capable of polymerizing to form a polymer having a high oxygen permeability", as used herein, refers to monomers, oligomers, macromers, and the like, and mixtures thereof, which are capable of polymerizing with like or unlike polymerizable materials to form a polymer which displays a relatively high rate of oxygen diffusion therethrough. For convenience of reference, these materials will be referred to herein as "oxyperm polymerizable materials" and the resultant polymers will be referred to herein as "oxyperm polymers".

The "oxygen transmissibility" of a lens, as used herein, is the rate at which oxygen will pass through a specific ophthalmic lens. Oxygen transmissibility, Dk/t , is conventionally expressed in units of barrers/mm, where t is the average thickness of the material [in units of mm] over the area being measured and "barrer" is defined as:

$$[(\text{cm}^3 \text{ oxygen})(\text{mm})/(\text{cm}^2)(\text{sec})(\text{mm Hg})] \times 10^{-9}$$

The "oxygen permeability", Dk , of a lens material does not depend on lens thickness. Oxygen permeability is the rate at which oxygen will pass through a material. Oxygen permeability is conventionally expressed in units of barrers, where "barrer" is defined as:

$$[(\text{cm}^3 \text{ oxygen})(\text{mm})/(\text{cm}^2)(\text{sec})(\text{mm Hg})] \times 10^{-10}$$

US 6,951,894 B1

5

These are the units commonly used in the art. Thus, in order to be consistent with the use in the art, the unit "barrer" will have the meanings as defined above. For example, a lens having a Dk of 90 barrers ("oxygen permeability barrers") and a thickness of 90 microns (0.090 mm) would have a Dk/t

A "polymerizable material which is capable of polymerizing to form a polymer having a high ion permeability", as used herein, refers to monomers, oligomers, macromers, and the like, and mixtures thereof, which are capable of polymerizing with like or unlike polymerizable materials to form a polymer which displays a relatively high rate of ion or water permeation therethrough. For convenience of reference, these materials will be referred to herein as "ionoperm polymerizable materials" and the resultant polymers will be referred to herein as "ionoperm polymers".

A "macromer", as used herein, refers to a polymerizable material which has a molecular weight of at least about 800 grams/mol. The term "macromer", as used herein, also encompasses oligomers.

A "monomer", as used herein refers to a polymerizable material which has a molecular weight of less than about 800 grams/mol.

A "phase", as used herein, refers to a region of substantially uniform composition which is a distinct and physically separate portion of a heterogeneous polymeric material. However, the term "phase" does not imply that the material described is a chemically pure substance, but merely that certain bulk properties differ significantly from the properties of another phase within the material. Thus, with respect to the polymeric components of a lens, an ionoperm phase refers to a region composed of essentially only ionoperm polymer (and water, when hydrated), while an oxyperm phase refers to a region composed of essentially only oxyperm polymer.

A "continuous phase" as used herein, refers to a region of substantially uniform composition which forms a continuous pathway from one surface of an article to another surface of an article.

"Co-continuous phases", as used herein, refers to at least two regions, each of substantially uniform composition which differs from the other, and each of which forms a continuous pathway from one surface of an article to another surface of an article. Thus, an ophthalmic lens having co-continuous phases of oxyperm polymer and ionoperm polymer will have two continuous pathways or sets of continuous pathways extending from the inner surface of the lens to the outer surface of the lens.

"Morphology", as used herein, refers to the structure and relationship of the phases of a material.

"Ophthalmically compatible" as used herein, refers to a material or surface of a material which may be in intimate contact with the ocular environment for an extended period of time without significantly damaging the ocular environment and without significant user discomfort. Thus, an ophthalmically compatible contact lens will not produce significant corneal swelling, will adequately move on the eye with blinking to promote adequate tear exchange, will not have substantial amounts of lipid adsorption, and will not cause substantial wearer discomfort during the prescribed period of wear.

"Ocular environment" as used herein, refers to ocular fluids (e.g., tear fluid) and ocular tissue (e.g., the cornea) which may come into intimate contact with a contact lens used for vision correction, drug delivery, wound healing, eye color modification, or other ophthalmic applications. "Hydrophilic" as used herein, describes a material or portion thereof which will more readily associate with water than with lipids.

6

A "hydrophilic surface" as used herein, refers to a surface which is more hydrophilic and lipophobic than the bulk or core material of an article. Thus, an ophthalmic lens having a hydrophilic surface describes a lens having a core material having a certain hydrophilicity surrounded, at least in part, by a surface which is more hydrophilic than the core.

The "outer surface" of a lens, as used herein, refers to the surface of the lens which faces away from the eye during wear. The outer surface, which is typically substantially convex, may also be referred to as the front curve of the lens. The "inner surface" of a lens, as used herein, refers to the surface of the lens which faces towards the eye during wear. The inner surface, which is typically substantially concave, may also be referred to as the base curve of the lens.

"TRIS" as used herein, refers to 3-methacryloxypropyltris(trimethylsiloxy) silane, which is represented by CAS No. 17096-07-0. The term "TRIS" also includes dimers of 3-methacryloxypropyltris(trimethylsiloxy) silane.

"Molecular weight" of a polymeric material (including monomeric or macromeric materials), as used herein, refers to the number-average molecular weight unless otherwise specifically noted or unless testing conditions indicate otherwise.

A. Oxyperm Polymerizable Materials

Oxyperm polymerizable materials include a wide range of materials which may be polymerized to form a polymer displaying a relatively high oxygen diffusion rate therethrough. In addition, these materials must be relatively ophthalmically compatible. These oxyperm polymerizable materials include, without limitation thereto, siloxane-containing macromers and monomers, fluorine-containing macromers and monomers, and carbon-carbon triple bond-containing macromers and monomers. The oxyperm macromer or monomer may also contain hydrophilic groups.

Preferred oxyperm polymers are those formed from a siloxane-containing macromer. Macromers having dialkyl siloxane groups, especially dimethyl siloxanes, are particularly preferred. These macromers are broadly referred to as poly(dimethyl siloxanes) (also, PDMS). The siloxane-containing macromer may also include hydrophilic groups. Examples of suitable siloxane-containing macromers include, without limitation thereto, the Materials A, B, C, and D as described herein.

The oxygen transmissibility (Dk/t) of the lens is preferably at least 70 barrers/mm, more preferably at least 75 barrers/mm, and most preferably at least 87 barrers/mm. The lens center thickness is typically more than about 30 microns, preferably about 30 to about 200 microns, more preferably about 40 to about 150 microns, even more preferably about 50 to about 120 microns, and most preferably about 60 to about 100 microns.

The oxygen transmissibility of the extended-wear lens from the outer surface to the inner surface must be sufficient to prevent any substantial corneal swelling during the period of extended wear. It is known that the cornea swells approximately 3% to 4% during overnight periods of sleep when the eyelids are closed, as a result of oxygen deprivation. It is also known that wearing a typical contact lens, such as ACUVUE (Johnson & Johnson), for a period of about 8 hours (overnight wear) causes corneal swelling of about 11%. However, a preferred extended-wear contact lens will produce, after wear of about 24 hours, including normal sleep periods, corneal swelling of less than about 8%, more preferably less than about 6%, and most preferably less than about 4%. A preferred extended-wear contact lens will produce, after wear of about 7 days, including normal sleep periods, corneal swelling of less than about 10%, more

US 6,951,894 B1

7

preferably less than about 7%, and most preferably less than about 5%. Thus, the extended-wear lens must have oxyperm polymer in an amount sufficient to produce oxygen diffusion pathways from the outer surface to the inner surface of the lens which are sufficient to yield the above properties relating to corneal swelling. Preferably, the extended-wear lens has a continuous phase of oxyperm polymer extending from the outer surface to the inner surface of the lens.

B. Ionoperm Polymerizable Materials

Ionoperm polymerizable materials include a wide range of materials which may be polymerized to form a polymer displaying a relatively high ion diffusion rate therethrough. In addition, these materials must be relatively ophthalmically compatible. These ionoperm polymerizable materials include, without limitation thereto, acrylates and methacrylates, such as 2-hydroxyethyl methacrylate, acrylamide, methacrylamide, and dimethylacrylamide; poly (alkylene glycols), such as poly(ethylene glycol); N-vinyl pyrrolidones such as N-vinyl-2-pyrrolidone; and the like and mixtures thereof. Other ionoperm materials are disclosed in the specific embodiments of Materials A-D, described below.

C. Weight Ratios

The ratios of oxyperm to ionoperm polymerizable materials may vary substantially, depending on the selected balance of oxygen permeability and ion permeability for the chosen end-use of the molded polymeric article. Preferably, the volumetric ratio of oxyperm to ionoperm material (including water) in the fully hydrated lens is about 40 to about 60 to about 60 to about 40. However, weight percentages, based on the total weight of the lens, will be defined because weight percentages are more conveniently utilized in lens fabrication. Preferably, the extended-wear contact lenses having substantially only ionoperm and oxyperm materials will have about 60 to about 85 weight percent oxyperm polymerizable material and about 15 to about 40 weight percent ionoperm polymerizable material in the prepolymerization mixture, based on total polymerizable material weight. More preferably, the prepolymerization mixture will contain about 70 to about 82 weight percent oxyperm polymerizable material and about 18 to about 30 weight percent ionoperm polymerizable material, based on total polymerizable material weight.

A wide variety of additional polymerizable materials may be included in the mixture prior to polymerization. Cross-linking agents, such as ethylene glycol dimethacrylate (EGDMA), may be added to improve structural integrity and mechanical strength. Antimicrobial polymerizable materials such as poly(quaternary ammonium) salts may be added to inhibit microbial growth on the lens material. Also, additional ionoperm monomers or macromers and oxyperm polymerizable materials may be added to adjust the oxygen permeability and ion permeability of the final molded article. An especially advantageous polymerizable material is TRIS, which may act both to increase oxygen permeability and to improve the modulus of elasticity.

A preferred prepolymerization mixture will include (a) about 30 to 60 weight percent oxyperm macromer, (b) about 20 to 40 weight percent ionoperm polymerizable material, and (c) about 1 to 35 weight percent TRIS, based on the total lens weight. More preferably, the amount of TRIS is about 10 to 33 weight percent, based on the total prepolymerization mixture weight.

In a preferred embodiment, the prepolymerization mixture includes less than about 5 weight percent cross-linking agent, based on the total prepolymerization mixture weight. More preferably, the prepolymerization mixture includes

8

less than about 2 weight percent cross-linking agent, based on the total prepolymerization mixture weight. Even more preferably, the prepolymerization mixture includes substantially no cross-linking agent. In a particularly preferred embodiment, the prepolymerization mixture includes no added cross-linking agent.

The previously described ranges for oxyperm polymerizable materials, ionoperm polymerizable materials, and TRIS are offered to enable the reader to better comprehend the invention. However, it should be noted that the specific weight or volume percentages of oxyperm and ionoperm polymerizable materials are not the most critical factors to consider in preparing a good extended-wear ophthalmic lens. More importantly, the lens must have sufficient ion permeability for good on-eye movement and sufficient oxygen permeability for good corneal health during the extended wear period.

D. Morphology

One requirement of the lens material is that the lens allow a high visible light transmission from the outer to inner surface of the lens. A lens morphology which includes large phase separated regions will reduce visible light transmission and cause substantial undesirable image distortion, thereby destroying the value of the lens as a vision correction device. Thus, the lens must have a morphology which allows at least about 80%, more preferably about 90%, visible light transmission and does not produce any significant undesirable image distortion.

In one preferred embodiment, the lens material has at least two phases: the phases including at least one oxyperm phase and at least one ionoperm phase. While there may be two distinct phases, it is believed that there may be a transition phase, or interphase, in which the material composition and the material properties are a blend of those of the oxyperm and ionoperm materials. Thus, there may exist a distinct oxyperm phase or plurality of distinct oxyperm phases, a distinct ionoperm phase or a plurality of distinct ionoperm phases, and an amphipathic phase mixture or blend of oxyperm and ionoperm phases. In one preferred embodiment, the glass transition temperature (T_g) of the oxyperm phase is less than about -115° Celsius.

The existence of separate oxyperm and ionoperm phases, rather than a complete blend of oxyperm and ionoperm phases, is believed to be advantageous in promoting the diffusion of oxygen and ions. Oxygen will diffuse predominantly through the oxyperm polymer, while the ionoperm polymer provides a higher barrier to oxygen diffusion. Similarly, ions will diffuse well through the ionoperm polymer, but the oxyperm polymer provides a higher resistance to ion diffusion. Thus, one homogeneous oxyperm/ionoperm phase will provide undesirable resistance to both oxygen and ion diffusion, while two separate oxyperm and ionoperm phases will provide low resistance pathways for transmission of both oxygen and ions or water. Thus, the ideal extended-wear lens has a pathway or series of pathways from the outer surface to the inner surface for transmission of oxygen therethrough, and an analogous continuous pathway or series of pathways for transmission of water or ions therethrough. In a particularly preferred embodiment, the lens has two co-continuous phases, one an oxyperm phase and the other an ionoperm phase, allowing for permeation of water or ions and oxygen between the front and base curves of the lens.

E. Bulk Water Content

The measurement of water content is difficult because the removal of adhered surface droplets, without affecting the

US 6,951,894 B1

9

bulk lens water content, is difficult. In addition, water may evaporate from the lens surface quickly, thereby lowering the water content from the equilibrium level. Accordingly, a discussion of the bulk water content of a lens warrants a discussion of the measurement technique used to determine the water content.

The preferred bulk water content of the hydrated lens will be a function of the lens material properties. The material properties are dependent on the prepolymerization macromers and monomers and polymerization conditions. Thus, the preferred water content for a lens including a fluorine-containing oxypm material may be different from that of a lens including a siloxane-containing oxypm material. Accordingly, while general ranges for bulk water content are offered for a better understanding of the invention, the invention is not generally limited to specific bulk water contents.

One method of measuring the water content of a lens formed in accordance with the present invention, referred to herein as the "Bulk Technique" is as follows. First the lens is thoroughly hydrated in a physiological saline solution, such that the water in the lens is in equilibrium with the surrounding water. Next the lens is gently blotted between two lint-free blotting cloths to remove surface moisture. The lens is quickly placed on an aluminum weighing pan and the first wet weight, W_1 , is measured. Next, the aluminum pan with lens is placed in an oven at 36° C. for a period of at least 24 hours. After heat treating, the pan with lens is removed, placed in a desiccator, and allowed to cool to room temperature (about 22° C.). The pan with lens is weighed again to determine the dry weight, W_d . The lens is re-equilibrated in physiological saline solution and a second wet weight W_2 is determined thereon. The wet weights (W_1 and W_2) are averaged to yield an average wet weight, W_w . The bulk water content is determined by the following equation:

$$\text{Percent water content} = (W_w - W_d) / W_w \times 100$$

A preferred lens bulk water content, determined by the "Bulk Technique" is less than about 32 weight percent. More preferably, the lens has a water content of about 10 to 30 weight percent, based on the total lens weight. A particularly preferred lens water content is about 15 to about 25 weight percent.

F. Ion and Water Permeability

Unexpectedly, it has been determined that the ion permeability through the lens correlates well with on-eye movement. As discussed earlier, it is known that on-eye movement of the lens is required to ensure good tear exchange, and ultimately, to ensure good corneal health. While the invention is not bound by theory presented herein, it may be useful to discuss some theory for a better understanding of ways to practice the invention.

It is theorized that water permeability is an exceptionally important feature for an extended-wear lens which includes oxypm polymers such as those disclosed herein. Siloxane-containing oxypm materials tend to adhere strongly to the eye, thereby stopping on-eye movement. The ability to pass water through the lens is believed to allow a siloxane-containing polymeric lens to move on the eye, where the movement occurs via forces exerted by water being squeezed out of the lens. The water permeability of the lens is also believed important in replenishing lens water content once pressure is removed. Further, the permeability of ions is believed to be directly proportional to the permeability of water. Thus, ion permeability is a predictor of on-eye movement.

However, regardless of whether the water permeability theory is a correct understanding of the actual on-eye

10

movement phenomenon, it has been unexpectedly found that above a certain threshold of ion permeability through a lens, from the inner surface of the lens to the outer, or vice versa, the lens will move on the eye, and below the threshold the lens will adhere to the eye. Thus, the present innovative extended-wear contact lenses provide a balance between the relatively high oxygen permeability (and associated high binding capacity) of oxypm materials with the low binding capacity (high on-eye movement) of ionopm materials. It is believed that this is accomplished by providing a plurality of continuous ion transmission pathways for ion and water movement through the lens.

It should be noted that ions may move through the lens via these ion pathways by a number of means. For example, ions may diffuse through the lens because of concentration differences from one surface to another. Ions may also be forced through the ion pathways by the mechanical action of blinking, with the concomitant compression forces on the lens essentially squeezing water out of the lens. In addition, the charge nature of the surfaces may provide an electromotive force which drives ion permeation through the lens. At times, one of these driving forces may be larger than the others, while at other times the relative magnitude may reverse. This discussion is presented to clarify that the invention is not restricted by the method or driving force by which ions move through the lens.

Neither the measurement of water permeability nor ion permeability through an ophthalmic lens is considered a routine matter of testing in the industry. Accordingly, a discussion of the preferred ion or water permeability ranges warrants a discussion of the measurement techniques used to determine the permeability.

The water permeability of a lens may be determined from the rate of water permeation through the lens, from one surface to another surface. The water permeability of a lens may be determined by positioning a lens between two reservoirs holding solutions having known, and different, initial concentrations of radiolabeled water (e.g., tritiated water), and then measuring concentration of radiolabeled water in the "receiving" reservoir (the reservoir towards which the net flow of radiolabeled water is positive) as a function of time.

The relative ion permeability of a lens may be determined from the rate of ion permeation through the lens, from one surface to another surface. The rate of ion permeation may be determined by positioning a lens between two reservoirs holding solutions having known, and different, initial ion concentrations, and then measuring conductivity in the "receiving" reservoir (the reservoir towards which the net flow of ions is positive) as a function of time. The concentration of ions, such as sodium, can be measured accurately using a pH meter and an ion-selective electrode. Ions are believed to be transmitted through a lens, from inner to outer surfaces and vice versa, primarily by the diffusion of ions through water pathways in the lens. Ion permeability through a lens is believed to be directly proportional to water permeability through a lens.

1. Ionoflux Measurement Technique

The following technique, referred to herein as the "Ionoflux Technique", is a preferred method for determining the ion permeability of a lens. This technique may be used to determine the likelihood of adequate on-eye movement.

The "Ionoflux Technique" involves the use of a conductometer (LF 2000/C, catalog no. 35 300105, Wissenschaftlich-Technische Werkstätten GmbH (WTW), Germany), an electrode equipped with a temperature sensor (LR 01/T, catalog no. 302 520, WTW), a donor chamber

US 6,951,894 B1

11

containing a salt solution, a receiving chamber containing about 60 ml of deionized water, a stir bar and a thermostat.

The donor chamber is specially designed for sealing a contact lens thereto, so that the donor solution does not pass around the lens (i.e., ions may only pass through the lens). The donor chamber is composed of a glass tube which is threaded at the end which is immersed in the receiving solution. The glass tube includes a centrally located hole of about 9 mm in diameter. A lid, which is threaded to mate with the glass tube, holds a lens-retaining member which includes a centrally located hole of about 8 mm in diameter. The lens-retaining member includes a male portion adapted to mate with and seal the edges of the inner (concave) surface of a lens and a female portion adapted to mate with and seal the edges of the outer (convex) surface of a lens.

The lens to be measured is placed in the lens-retaining member, between the male and female portions. The male and female portions include flexible sealing rings which are positioned between the lens and the respective male or female portion. After positioning the lens in the lens-retaining member, the lens-retaining member is placed in the threaded lid. The lid is screwed onto the glass tube to define the donor chamber. The donor chamber is filled with 16 ml of 0.1 molar NaCl solution. The receiving chamber is filled with 60 ml of deionized water. The leads of the conductivity meter are immersed in the deionized water of the receiving chamber and a stir bar is added to the receiving chamber. The receiving chamber is placed in a thermostat and the temperature is held at about 35° C. Finally, the donor chamber is immersed in the receiving chamber.

Measurements of conductivity are taken every 20 minutes for about three hours, starting 10 minutes after immersion of the donor chamber into the receiving chamber. The Ionoflux Diffusion Coefficient, D, is determined by applying Fick's law as follows:

$$D = -n' / (A \times dc / dx)$$

where

n' = rate of ion transport [mol/min]

A = area of lens exposed [mm²]

D = Ionoflux Diffusion Coefficient [mm²/min]

dc = concentration difference [mol/L]

dx = thickness of lens [mm]

An Ionoflux Diffusion Coefficient of greater than about 1.5×10^{-6} mm²/min is preferred for achieving sufficient on-eye movement. More preferably, the Ionoflux Diffusion Coefficient is greater than about 2.6×10^{-6} mm²/min, while most preferably, the Ionoflux Diffusion Coefficient is greater than about 6.4×10^{-6} mm²/min. It must be emphasized that the Ionoflux Diffusion Coefficient correlates with ion permeability through the lens, and thereby is a predictor of on-eye movement.

2. Ionoton Measurement Technique

The following technique, referred to herein as the "Ionoton Technique", is another preferred method for determining the relative ion permeability of a lens. The technique is based on measurement of the diffusion of sodium chloride through a lens.

The "Ionoton Technique" involves the use of a pH meter (Beckman, VWR catalog no. BK123142), a VSC-1 Diffusion Cell Drive Console (Crown-Bio, Somerville, N.J.), a DCB-100B Diffusion Cell (Crown-Bio), and a 6 cm sodium ion-specific electrode (Microelectronics, Londonderry, N.H., catalog no. MI-414P). The technique is not limited to the aforementioned instruments or materials; equivalent instruments or materials may be used.

12

First, a contact lens is mounted onto an orifice of the DCB-100B cell chamber, the donor chamber. Next, the connecting cell chamber (receptor chamber) is placed against the cell chamber containing the contact lens and tightly clamped on the clamp holder supplied with the VSC-1 Drive Console. Then, a phosphate-buffered saline (PBS, Mediatech catalog no. 21-031-LV) is placed into the receptor side of the cell chamber. Stir bars are added to each cell chamber. The 6 cm electrode is placed into the PBS saline receptor side. After the electrode has equilibrated in the PBS saline, the pH meter is placed in the mV function to establish the 0 mV point. PBS which has been saturated with sodium chloride is added to the donor chamber.

The millivolt signal is recorded at 5, 10, 15, 30, 60, 120, and 180 minute intervals. The millivolt signal is converted to a sodium ion concentration by a standard curve of sodium ion concentration vs. millivolt signal. The Ionoton Ion Permeability Coefficient, P, is then determined in accordance with the following equation:

$$\ln(1 - 2C(t)/C(0)) = 2APt/Vd$$

where:

$C(t)$ = concentration of sodium ions at time t in the receiving cell

$C(0)$ = initial concentration of sodium ions in donor cell

A = membrane area, i.e., lens area exposed to cells

V = volume of cell compartment (3.0 ml)

d = average lens thickness in the area exposed

P = permeability coefficient

The average thickness of the lens in the exposed test area may be determined by averaging a number of readings, e.g., 10 readings, with a low-pressure thickness-measuring instrument, such as a Mitotoya micrometer VL-50, or equivalents thereof. The Ionoton Ion Permeability Coefficient, P, having units of cm²/second, may be determined from the slope of a plot of time (t) v. $\ln(1 - 2C(t)/C(0)) \times (-2At/Vd)$.

An Ionoton Ion Permeability Coefficient, P, of greater than about 0.2×10^{-6} cm²/sec is preferred, while greater than about 0.3×10^{-6} cm²/sec is more preferred and greater than about 0.4×10^{-6} cm²/sec is most preferred. It must be emphasized that the Ionoton Ion Permeability Coefficient correlates with ion permeability through the lens, and thereby is a predictor of on-eye movement.

3. Hydrodell Water Permeability Technique

The following technique, referred to herein as the "Hydrodell Technique" is a preferred method for determining the water permeability of a lens. This technique may be used to determine the likelihood of adequate on-eye movement.

The Hydrodell Technique involves the measurement of the rate of transfer of the radiolabeled solutes THO (³H—HO or tritiated water) and ¹⁴C-glucose across the contact lens using a two-chamber apparatus. ¹⁴C-glucose is used in this measurement technique to reveal any leak in the system during testing. The lens is mounted between chambers, which are stirred at a controllable rate. Chamber I contains a solution with a high concentration of labeled solute. Chamber II, the "receiving chamber" contains an identical solution but without the labeled solute. Samples of the solution in chambers I and II are taken at intervals over the test period. The radioactivity in the samples is measured. The permeability of the lens is calculated from the measured radioactivity, the sample times, the chamber volumes and the lens area exposed to the solutions. A more detailed description of the Hydrodell Technique follows.

US 6,951,894 B1

13

a. Solution Preparation

Dulbecco's phosphate buffered saline (DPBS) is prepared by first dissolving, sequentially, about 160 g sodium chloride (NaCl), about 4 grams potassium chloride (KCl), about 23 grams disodium hydrogen orthophosphate (Na_2HPO_4), about 4 grams potassium dihydrogen orthophosphate (KH_2PO_4), and about 10 grams sodium azide in a liter of reverse-osmosis (MilliQ) water. Then, the pH is adjusted to about 7.3 by adding appropriate amounts of HCl. Finally, the buffer solution is diluted to 1:20 (50 ml buffer solution with 950 ml reverse-osmosis water), and allowed to degas either in a screw-capped container overnight or under vacuum.

A Cold Glucose buffer solution is prepared by adding about 0.1 grams D-glucose to one liter of DPBS, followed by sterilization via filtration through a 0.2 μl millipore filter and storage at 4° C. until use.

The Chamber I solution is prepared by adding about 6 μl THO (TR53, 1.0 mCi/ml activity, available from Amersham Australia, located in North Ryde NSW Australia) and about 16 μl ^{14}C -glucose THO (in ethanol, available from Amersham Australia) to about 12 ml of the Cold Glucose buffer solution. Preferably, this solution is used within about 24 hours of preparation. The Chamber II solution is DPBS.

b. Apparatus Preparation

The chambers have a volume sufficient to hold about 12 ml of solution during testing. While the exact shape of the chambers is not critical, both chambers have rectangular cross-sections for ease of construction. The chambers may be made from a variety of water-proof rigid materials, preferably clear (e.g., acrylic plates, FX Plastics, Marrickville NSW Australia) so that samples may be observed during testing. Each chamber has a circular aperture of about 7 mm diameter appropriate for mounting a lens between the chambers for contact with solutions held within the chambers. Some affixing means, such as a set of mounting bolts, are necessary to securely affix one chamber to the other with the lens mounted in between.

A test contact lens is mounted symmetrically over the aperture of Chamber II. Folds and wrinkles are manually removed from the lens. Chamber I is positioned adjacent the aperture and mounted lens of Chamber II, and the chambers are secured to one another using mounting bolts.

About 12 ml (V_2) of DPBS is placed in Chamber II. About 12 ml of the Chamber I labeled solution is placed in Chamber I, at which point time $t=0$ is established. Stirrers are added to both chambers and the stirrer speed is set at about 1200 rpm.

c. Sampling

Sampling generally starts at time $t_0=5$ minutes. The final sample time, t_p , is usually at about 50 minutes for high water content lenses and about 120 minutes for low water content lenses, although these times are not critical.

At time $t_0=5$ minutes, two samples of about 0.2 ml volume are pipetted from Chamber I, and two 0.2 ml aliquots of DPBS are added to Chamber I to restore the volume. These samples are placed into plastic counting tubes with about 4 ml Ultima Gold™ cocktail (available from Packard Instrument Co., Meriden, Conn.) and about 0.9 ml DPBS.

Also at time t_0 , one sample of about 1.0 ml volume is pipetted from Chamber II and one 1.0 ml aliquot of DPBS is added to Chamber II to restore the volume. The sample is placed into a plastic counting tube with about 4 ml Ultima Gold™ cocktail.

At intermediate times between t_0 and t_p (e.g., every 10 minutes), one sample of about 1.0 ml volume is pipetted from Chamber II and one 1.0 ml aliquot of DPBS is added to Chamber II to restore the volume. Each sample is placed into a plastic counting tube with about 4 ml Ultima Gold™ cocktail.

14

At time t_p , two samples of about 0.2 ml volume are pipetted from Chamber I. These samples are placed into plastic counting tubes with about 4 ml Ultima Gold™ cocktail and about 0.9 ml DPBS.

Also at time t_p , two samples of about 1.0 ml volume are pipetted from Chamber II. These samples are placed into plastic counting tubes with about 4 ml Ultima Gold™ cocktail.

d. Measurements

The activity of the samples are measured by liquid scintillation counting, or other appropriate technique. Liquid scintillation counting may be advantageously accomplished by using protocol number 6 for $^3\text{H}/^{14}\text{C}$ on a Tri-Carb Liquid Scintillation Analyzer (1900TR, available from Packard Instrument Co.).

Three standards containing about 10^4 to 10^5 cpm THO in reversed-osmosis (MilliQ) water are prepared. Three standards containing about 10^4 to 10^5 cpm ^{14}C glucose in reversed-osmosis (MilliQ) water are also prepared. A blank containing MilliQ water is prepared.

The scintillation analyzer settings are LLA=0 KeV and ULA=12 KeV for ^3H ("1") in channel 1 and LLB=12 KeV and ULB=156 KeV for ^{14}C ("2") in channel 2. The standards and blank are counted three times during each counting of samples, and the counts are averaged. The following denote the relevant measured sample activities:

b_1 =measured activity of blank sample in channel 1

b_2 =measured activity of blank sample in channel 2

S'_{11} =measured activity of standard ^3H sample in channel 1

S'_{12} =measured activity of standard ^{14}C sample in channel 2

S'_{21} =measured activity of standard ^3H sample in channel 1

S'_{22} =measured activity of standard ^{14}C sample in channel 2

y_1 =measured activity of test sample (both ^3H and ^{14}C) in channel 1

Y_2 =measured activity of test sample (both ^3H and ^{14}C) in channel 2

e. Water Permeability Calculation

In order to calculate the actual activity of a sample, the measured activities of the isotopes, ^3H and ^{14}C , must first be corrected to remove the cross-contamination error due to the presence of both isotopes in one sample. Without explaining the mathematical derivations, the following stepwise procedure is offered as an example of one method of determining water permeability from the above measurements:

(1) Calculate S_{11} , S_{12} , S_{21} , and S_{22} , from the following equations:

$$S_{11}=S'_{11}-b_1$$

$$S_{12}=S'_{12}-b_1$$

$$S_{21}=S'_{21}-b_2$$

$$S_{22}=S'_{22}-b_2$$

(2) Calculate a_{12} and a_{21} from the following equations:

$$a_{12}=S_{12}/S_{22}$$

$$a_{21}=S_{21}/S_{11}$$